CONSENSUS ORTHOPEDICS INC. CONSENSUS® REVISION KNEE SYSTEM Instructions for Use (IFU)

IMPORTANT INFORMATION FOR SURGEON: PLEASE READ PRIOR TO IMPLANTING THIS DEVICE IN A CLINICAL SETTING. THE SURGEON SHOULD BE FAMILIAR WITH THE SURGICAL TECHNIQUE.

DESCRIPTION The CONSENSUS® REVISION KNEE SYSTEM (RKS) contains left and right configurations of a femoral component and fixed-bearing tibial baseplate... PATIENT COUNSELING INFORMATION Complications and/or failure of revision knee prostheses are more likely to occur in patients with unrealistic functional expectations...

HOW PRODUCT IS SUPPLIED Each component of the CONSENSUS® RKS is supplied in individually sterilized boxes or packages, in a wide range of sizes. Please refer to the current price list, surgical technique, or catalog for the catalog number and sizes available.

INFORMATION FOR USE The recommended trial components are used for size determination, preparation, evaluation, trial reduction, and range of motion evaluation. The use of trials preserves the integrity of implants and sterile packaging. Trials are not intended for implantation.

INFORMATION FOR USE The indications for the use of the CONSENSUS® RKS must be carefully considered with respect to the patient's entire evaluation and alternative procedures. Patient selection is dependent on age, general health, available bone stock and quality, and any prior surgery or anticipated further surgery.

ASSEMBLY OF COMPONENTS The CONSENSUS® RKS tibial component consists of two parts to be assembled at the time of surgery, the CONSENSUS® RKS tibial baseplate and either the congruent, ultracongruent, or PS configurations of the CONSENSUS® RKS femoral component.

Table with 6 columns: Component, Size 1, 2, 3, 4, 5, 6. Rows include Femoral Component, Femoral Posterior Augment (5mm), Femoral Distal Augment (5mm), Knee Tibial Insert, Knee All Poly Tibia, Knee Patella, Knee Patella Oval, Knee Patella Oval Metal Backed, and MKB Tibial Insert.

Table with 7 columns: Component and Size Compatibility with CONSENSUS® RKS Tibial Baseplate, Tibial Baseplate, Tibial Augment (5mm), Tibial Augment (10mm), Tibial Augment (15mm), Tibial Augment (20mm), Tibial Augment (25mm).

Table 1: Components compatible with the CONSENSUS® RKS femoral and tibial components. COI recommends to follow the size compatibility table for optimal surgical outcome.

- INDICATORS FOR USE A. Primary intervention of rheumatoid arthritis, osteoarthritis, traumatic arthritis, polyarthritis, collagen disorders, and/or avascular necrosis of the femoral condyle. B. Post-traumatic loss of joint configuration (particularly when there is patellofemoral erosion, dysfunction, or prior patellectomy).

INTENDED PERFORMANCE All components of the CONSENSUS® RKS are intended to perform in a safe and effective manner in restoring knee function within the intended use of the product.

- CONTRAINDICATIONS A. Any active or suspected latent infection in or about the knee joint. B. Mental or neuromuscular disorders which would create unacceptable risk of prosthesis instability or complications in post-operative care.

- PRECAUTIONS A. Before any implant is used the surgeon should be completely familiar with all aspects of the surgical procedure and the limitations of the device. B. Instruct patients on the limitations of the prosthesis and how to modify their activities accordingly.

- WARNINGS A. All CONSENSUS® RKS components are supplied sterile. The physician's insert and double Tyvek packages are subject to ethylene oxide sterilization. Do not use any component if the package has been breached.

- ADVERSE EFFECTS A. All prosthetic replacements have the potential for adverse effects, including infection, loosening, fracture, wear, and/or fracture of the components, disassembly, or positional changes of the components. B. Systemic reactions to component materials could occur, and should be ruled out preoperatively.

- or tissue inflammation of unknown etiology. E. Reoperation may be necessary to correct adverse effects. F. On rare occasions, complications may require arthrodesis or amputation of the limb.

CONSENSUS® REVISION KNEE SYSTEM (RKS) contains left and right configurations of a femoral component and fixed-bearing tibial baseplate. In addition to accessory components including stems, a taper plug, a baseplate plug, femoral and tibial augments, screws, and pegs, it is a non-constrained or semi-constrained design (depending on the insert selected), in which both medial and lateral collateral ligaments must be intact. All RKS stems are intended FOR PRESS-FIT NON-CEMENTED USE.

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Table with 6 columns: Component, Size 1, 2, 3, 4, 5, 6. Rows include Femoral Component, Femoral Posterior Augment (5mm), Femoral Distal Augment (5mm), Tibial Insert, Tibial Augment (5mm), Tibial Augment (10mm), Tibial Augment (15mm), Tibial Augment (20mm), Tibial Augment (25mm).

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